# DESS LAB

## IFU

CUSTOM ABUTMENT FOR DENTAL IMPLANTS

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DI-MUL-o8-DSL

03-04-2018

# Custom Abutment for dental implants

## INSTRUCTIONS FOR USE

ENGLISH (INTERNATIONAL)

### CAUTION – PLEASE READ CAREFULLY

DESS LAB® Custom Abutments should only be used by dental specialists with experience in maxillary implantology and other specialities, such as dental diagnosis, planning, dental surgery and CAD prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. DESS LAB® Custom Abutments and DESS® screws designed for single-use only must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health like patient tissue infection. All DESS LAB® Custom Abutments must be dry fitted before use to check the correct fitting. The clinician will be responsible for correct application of those restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor. **U.S. Federal law restricts these devices to sale by or on the order of a licensed dentist / physician.** We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must a screw check. If the screws have suffered unusual wear, the complete integrity of the implant-abutment should be checked. New screws must be used for any revision, correction or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty. **During any intraoral use and manipulation all DESS LAB® and DESS® products must be secured to prevent aspiration due to its shape and size.** The responsibility for the prosthesis and its construction that will be placed above DESS LAB® Custom Abutment shall lie with the person who has designed the *st/* file and who will complete the restoration because in any moment there is no way of interfering in abutment design on the part of DESS LAB® employees and the file remains unchanged, neither design validation is proceed. Please follow "DESS LAB DESIGN LIMIT PARAMETERS" brochure available at [www.dess-lab.com/downloads/1/download](http://www.dess-lab.com/downloads/1/download)

### INDICATIONS

Custom Implant Abutments are used for prosthetic restorations of dental implants. They are patient-specific product designed by Dental Laboratory using CAD Software. To use this product in digital dentistry DESS LAB® library and scan abutment are needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-lab.com/> After registration we will send you a link to download the files. Some dental CAD/CAM systems require a double scan with diagnostical wax-up. If you need additional assistance with installation and/or proper use of the files and with use of piece itself, please contact your DESS® local distributor or find necessary information in corresponding instructions. Designing the abutment follow the recommendations given by your ceramic provider and CAD technique design rules. Please follow "DESS LAB DESIGN LIMIT PARAMETERS" brochure available at [www.dess-lab.com/downloads/1/download](http://www.dess-lab.com/downloads/1/download)

### CONTRAINDICATIONS

All materials used are biocompatible; however, some patients may present allergies or hypersensitivity to any of the materials and its components (specified in the table). Do not use for restorations with cantilever on a single implant, bruxism or insufficient space.

### COMPATIBILITY INFORMATION

All DESS LAB® Custom Abutments are available for a variety of connections. For compatibility with dental implants and implant analogues please visit our website <http://www.dess-lab.com> or contact your local distributor. To find necessary information about ANKYLOS® C/X system, please read the last paragraph of these instructions.

### STERILIZATION AND RE-USE

All products are supplied NON STERILE. For sterilisation, we recommend autoclaving the product at 121°C for 30 minutes, drying time 30 minutes (in accordance with standard UNE-EN ISO 17665-1:2007). They are marked for "Single use only" because it is a custom abutment designed for specific patient. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use custom device will not be accepted.

ARTICLE	GROUP DESS LAB REF. N°	MATERIAL	STERILIZATION	RE-USE
Co-Cr CUSTOM ABUTMENT	51.0XX-C-XX	Cobalt-Chromium-Molybdenum alloy	Autoclave before use on the patient	Single use <input checked="" type="checkbox"/>
TITANIUM CUSTOM ABUTMENT	51.0XX-T-XX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use <input checked="" type="checkbox"/>
TITANIUM CUSTOM ABUTMENT WITH MetAlive® COATING	51.0XX-M-XX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use <input checked="" type="checkbox"/>
DESS SCREW	19.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use <input checked="" type="checkbox"/>

### APPLICATIONS AND DIRECTIONS FOR USE – CoCr CUSTOM ABUTMENT

COMPOSITION (percentage by mass)	TECHNICAL DATA
Co 65.4 % Cr 27.75 % Mo 5.06 % Additional elements less than 1 % (Fe, Mn, Si, Ni, N, C, Ti, P, S)	Yield strength <sup>1</sup> Rp 0.2 > 827 MPa Tensile strength <sup>1</sup> Rm > 1172 MPa Hardness <sup>1</sup> HV10 459 Fracture elongation <sup>1</sup> A5 21 % Solidus <sup>1</sup> 1369 °C Liquidus <sup>1</sup> 1415 °C CTE <sup>1</sup> (25 – 500 °C) 14.1 x 10 <sup>-6</sup> K <sup>-1</sup> Color white

The Co-Cr Custom Abutment is used for prosthetic restorations prepared by dental technicians in the dental laboratory. It is an implant patient-specific device that consists of a mechanized abutment for direct ceramic firing by ceramist. For any laboratory process the abutment must be carefully hand-tightened into the lab analogue with a lab screw. The existing abutment screw should be used only for the final integration of the restoration after being sterilized.

Before ceramic layering remember to clean the structure with vapour and degrease the surface. If sandblasting is needed or neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. The veneering should have uniform thickness according to ceramic instructions to prevent cracks and chips.

Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend a radiography to verify the correct settlement after tightening. Screw should be fixed with adequate torque according to implant connection and size. Sealing of access hole with Teflon tape or gutta-percha and composite is needed. After placement remember to check static and dynamic occlusion using conventional protocol.

<sup>1</sup> Small variations between batches of sourced raw material are possible.

### APPLICATIONS AND DIRECTIONS FOR USE – TITANIUM CUSTOM ABUTMENT

The Titanium Custom Abutment is used for prosthetic restorations prepared by dental technicians in the dental laboratory. It is an implant patient-specific device that consists of a mechanized abutment forming the core of prosthetic crowns and bridges and used for their support. For any laboratory process the abutment must be carefully hand-tightened into the lab analogue with a lab screw. The existing abutment screw should be used only for the final integration of the restoration after being sterilized. If sandblasting is needed or neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts and that they are in the correct place and position for all secondary parts. We strongly recommend a radiography to verify the correct settlement after tightening. To improve cement adhesion we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (dual- or self-cured) intra-orally following recommendations of cement manufacturer. Before cementing the final restoration, check the torque according to implant connection and size. When cementing use definitive luting cement or use temporary cement after sealing of access hole with Teflon tape or gutta-percha. Remove cement excess from edge of the crown to avoid peri-implantitis which can lead to implant loss. After cementation remember to check static and dynamic occlusion using conventional protocol.

### APPLICATIONS AND DIRECTIONS FOR USE – METALIVE® COATING

Some DESS LAB® Titanium Custom Abutments are covered with MetAlive® coating which enables protein-mediated attachment of the gingival tissue to the abutment surface. However, this does not change the general indications for use of the modified abutment.

**Warning:** Do not treat the gingival contact surface with any bonding agents, cements, cleaning agents or disinfectants. Do not clean the surface using abrasive technologies. Plaque removal should be done using soft scaling instruments. If sandblasted in laboratory, the gingival contact area must be protected with wax. Before sterilization or for cleaning, water vapour can be used as well as acetone and absolute alcohol in ultrasonic cleaner for 5 min. After cleaning or sterilizing do not touch surface with bare fingers and keep the surface hygienic. Please avoid disinfectants such as iodine, Hydrogen peroxide, chlorhexidine, polyvidone, etc. Avoid pointing strong laser or IR-radiation to modified surface. Do not remove MetAlive® surface material from gingival contact area with dental burr or by polishing instruments such as diamond rubber. Be sure that the surface is not in contact with hard sharp edges or points. **Processing:** MetAlive® abutments can be sterilized in autoclave at 121° C. After sterilization have the gingival tissue wound refreshed to enhance tissue repair and attachment, also after replacement of healing abutment with the permanent MetAlive® abutment. Ensure maximal initial contact with gingival tissue. Ensure minimal movement of gingival tissue during healing. Patient should use soft toothbrushes and non-abrasive toothpaste. For cleaning the interdental spaces dental floss or interdental brushes can be used.

**APPLICATIONS AND DIRECTIONS FOR USE – DESS SCREW**


The screw is for fixing prosthetic restorations and auxiliary abutments over implant or analogue. Make sure to secure the parts with corresponding screw and observe specified torque value placed on the label or ask your local DESS® distributor. You can download torques brochure from <http://www.dess-abutments.com/downloads>. For best results, the following conditions must be meticulously met:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The driver should be placed in the longitudinal axle of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every checks thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case.
- If you have any doubts about the fit of the screw, consult the local distributor.
- **Position the patient to avoid aspiration in case the screw falls during screwing/unscrewing.**

**INTERNAL ANK SPECIAL REQUIREMENTS**

The thickness of the anti-rotation elements is reduced due to the widened screw channel in order to accommodate the use of a normal screw. To prevent bending of those parts under applied pressure, please ensure that the screw is threaded through the abutment whilst restoration takes place and before placing it in the analog or implant. If this precaution is taken, the anti-rotation elements remain supported and safe from shear and compression forces. If, for any reason the screw needs to be removed, place it back again in its position before applying any loads on the abutment.

**EXPLANATION OF THE SYMBOLS THAT APPEAR ON THE LABELS**

 For single use only	 Serial No.	 Reference
 NON Sterile product	 Read instructions	 Patient identification
 Manufacturer information	Manufacturer: Terrats Medical SL Avda. La Ferrería, 62- 08110 – Montcada i Reixac -Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 <a href="mailto:info@dess-lab.com">info@dess-lab.com</a>	

All DESS LAB® Custom Abutments are manufactured according to ISO 9001:2008 and ISO 13485:2012.

3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA  
 ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE®, FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants  
 BioHorizons® is a registered trademark of BioHorizons, INC.  
 Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and Multiunit® are registered trademarks of Nobel Biocare AB, Sweden  
 ITI®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland  
 MIS® is a registered trademark of MIS Implant Technologies Ltd.  
 Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA  
 Dyna® is a registered trademark of Dyna Dental Engineering BV  
 Locator® is a registered trademark of Zest Anchors INC.

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. DESS LAB® Custom Abutments may not be available in your country. These instructions replace all previous editions. For any detailed information about DESS LAB® products please contact your local distributor.

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ITI®, Straumann® et synOcta® sont des marques déposées de Straumann Holding AG, Suisse.
MIS® est une marque déposée de MIS Implant Technologies Ltd..